

# PierenFemme

## SIMPLIEST ELECTRO-STIMULATOR FOR THE TREATMENT OF INCONTINENCE

Millions of people are suffering from incontinence. Treatment with PierenFemme significantly improves the strength of the pelvic floor muscles and sphincter muscles.

PierenFemme is a single channel stimulator which includes 4 pre-set programs for each kind of patient:

- Urge incontinence
- Stress incontinence
- Mixed incontinence 1
- Mixed incontinence 2

PierenFemme can be used with vaginal or anal probes but is also effective with self-adhesive electrodes.



### SUPPLIED WITH THE FOLLOWING ACCESSORIES

462126	1	PierenFemme
106351	1	Cable type 5.15
103025	1	Vaginal probe V2B
450780-11	2	AAA batteries
451600-0336	1	Operating instructions

### TECHNICAL DATA

Item no.	462126	Output voltage	0 – 60 Vpp (on 500 $\Omega$ real)
Frequency			10 Hz to 50 Hz
Pulse width			250 $\mu$ s
Power supply			2 x 1.5V AAA batteries
Modulation type			Amplitude modulation with "impulse envelope curve"
Dimensions			120 x 34 x 20 mm
Weight			approx. 40 g without batteries